

Package ‘HiTMaP’

March 17, 2021

Type Package

Title A collection of tools for imaging MS data processing

Version 1.0.0

Date

Depends R (>= 4.0.2)

Author George GUO [aut, cre]

Maintainer George GUO <George.GUO@auckland.ac.nz>

Description

License GPL-3

LazyLoad yes

NeedsCompilation no

Encoding UTF-8

BugReports

Imports pryr,colortools,pacman,purrr,lattice,pls,
RColorBrewer,magick,ggplot2,reticulate,
dplyr,stringr,data.table,doParallel,
iterators,foreach,protViz,cleaver,
MALDIquant,Biostrings,XVector,IRanges,Cardinal,
ProtGenerics,S4Vectors,stats4,egg,
EBImage,BiocParallel,BiocGenerics,parallel,
stats,graphics,grDevices,utils,
datasets,methods,Rdisop,Rcpp,
pacman,devtools,fs,processx,S4Vectors,
spdep,zoo,KEGGREST,OrgMassSpecR,
enviPat,rgl,multtest,XML,
BiocManager,shinyFiles,shiny,
rcdk,colorspace,OneR,
gridExtra,grid,gtable,ChemmineR

VignetteBuilder knitr

Suggests knitr,
rmarkdown,
devtools,
plotly,
UniProt.ws,
bitops

RoxygenNote 7.1.1

biocViews Software, Infrastructure, Proteomics, Lipidomics,
MassSpectrometry, ImagingMassSpectrometry, ImmunoOncology,
Normalization, Clustering, Classification, Regression

R topics documented:

Cleavage_rules_fun	2
cluster_image_grid	3
get_atoms	3
grid.ftable	4
HiTMaP_GUI	4
imaging_identification	5
imaging_Spatial_Quant	9
merge_atoms	11
Meta_feature_list_fun	12
Peptide_modification	13
Protein_feature_list_fun	13
Index	15

Cleavage_rules_fun	<i>Cleavage_rules_fun</i>
--------------------	---------------------------

Description

This is a function will return a list of pre-defined enzyme digestion specificity.

Usage

Cleavage_rules_fun()

Value

a table of Cleavage rules

Examples

Cleavage_rules_fun()

cluster_image_grid	<i>cluster_image_grid</i>
--------------------	---------------------------

Description

This function renders the clustered images for maldi imaging data set

Usage

```
cluster_image_grid(  
  clusterID,  
  SMPLIST,  
  imdata,  
  ClusterID_colname = "Protein",  
  componentID_colname = "Peptide",  
  Component_plot_threshold = 2,  
  Component_plot_coloure = c("mono", "as.cluster"),  
  smooth.image = "gaussian",  
  contrast.enhance = "suppression",  
  colorpallet = "Set1",  
  plot_layout = c("line", "grid"),  
  export_Header_table = F,  
  export_footer_table = F,  
  plot_style = c("fleximaging", "ClusterOnly", "rainbow"),  
  protein_coverage = F,  
  footer_style = "Length",  
  output_png_width_limit = 1980,  
  attach_summary_cluster = T,  
  cluster_color_scale = c("blackwhite", "fleximaging"),  
  remove_cluster_from_grid = T,  
  img_brightness = 100,  
  ppm = 20,  
  list_of_protein_sequence,  
  workdir = getwd(),  
  pixel_size_um = 50,  
  Score_thres = NULL  
)
```

get_atoms	<i>get_atoms</i>
-----------	------------------

Description

This is a function that prepare the atoms list for candidate molecules.

Usage

```
get_atoms(Symbol)
```

Arguments

Symbol the Symbol of the molecule

Value

a atoms list

Examples

```
get_atoms(Symbol="C7H13O8P")
```

grid.ftable

grid.ftable

Description

This is a function will plot table.

Usage

```
grid.ftable(d, padding = unit(4, "mm"), ...)
```

Value

a table image

Examples

```
grid.ftable()
```

HiTMaP_GUI

HiTMaP shiny-based graphical user interface (GUI)

Description

HiTMaP shiny-based graphical user interface (GUI)

Usage

```
HiTMaP_GUI(wd = "", port = 3838)
```

Examples

```
## Start the GUI
## Not run:
HiTMaP_GUI()

## End(Not run)
```

`imaging_identification`*imaging_identification*

Description

This is a peptide mass fingerprint search function for maldi imaging data analysis

Usage

```
imaging_identification(  
  datafile,  
  projectfolder = NULL,  
  threshold = 0.005,  
  ppm = 5,  
  mode = c("Proteomics", "Metabolomics"),  
  Digestion_site = "trypsin",  
  missedCleavages = 0:1,  
  Fastadatabase = "uniprot-bovin.fasta",  
  adducts = c("M+H"),  
  Modifications = list(fixed = NULL, fixmod_position = NULL, variable = NULL,  
    varmod_position = NULL),  
  Substitute_AA = NULL,  
  Decoy_search = TRUE,  
  Decoy_adducts = c("M+ACN+H", "M+IsoProp+H", "M+DMSO+H", "M+Co", "M+Ag", "M+Cu",  
    "M+He", "M+Ne", "M+Ar", "M+Kr", "M+Xe", "M+Rn"),  
  Decoy_mode = "isotope",  
  mzrange = "auto-detect",  
  Database_stats = F,  
  adjust_score = FALSE,  
  IMS_analysis = TRUE,  
  PMFsearch = IMS_analysis,  
  Load_candidatelist = IMS_analysis || plot_cluster_image_grid,  
  Bypass_generate_spectrum = FALSE,  
  peptide_ID_filter = 2,  
  Protein_feature_summary = TRUE,  
  Peptide_feature_summary = TRUE,  
  plot_ion_image = FALSE,  
  parallel = detectCores(),  
  spectra_segments_per_file = 4,  
  Segmentation = c("spatialKMeans", "spatialShrunkenCentroids", "Virtual_segmentation",  
    "none", "def_file"),  
  Segmentation_def = "Segmentation_def.csv",  
  Segmentation_ncomp = "auto-detect",  
  Segmentation_variance_coverage = 0.8,  
  preprocess = list(force_preprocess = FALSE, use_preprocessRDS = TRUE, smoothSignal =  
    list(method = "disable"), reduceBaseline = list(method = "locmin"), peakPick =  
    list(method = "adaptive"), peakAlign = list(tolerance = ppm, units = "ppm"),  
    normalize = list(method = c("rms", "tic", "reference")[1], mz = 1)),  
  Smooth_range = 1,  
  Virtual_segmentation_rankfile = NULL,
```

```

Rotate_IMG = NULL,
Region_feature_summary = FALSE,
Spectrum_validate = TRUE,
output_candidatelist = TRUE,
use_previous_candidates = FALSE,
score_method = "SQ RTP",
plot_cluster_image_grid = FALSE,
plot_cluster_image_maxretry = 2,
plot_cluster_image_overwrite = F,
smooth.image = "gaussian",
componentID_colname = "Peptide",
ClusterID_colname = "Protein",
Protein_desc_of_interest = ".",
Protein_desc_of_exclusion = NULL,
plot_unique_component = TRUE,
FDR_cutoff = 0.05,
use_top_rank = NULL,
plot_matching_score = F,
Component_plot_coloure = "mono",
cluster_color_scale = "blackwhite",
plot_layout = "line",
export_Header_table = T,
export_footer_table = T,
attach_summary_cluster = T,
remove_cluster_from_grid = attach_summary_cluster,
pixel_size_um = 50,
img_brightness = 100,
Thread = NULL,
cluster_rds_path = NULL,
remove_score_outlier = F,
Plot_score_IQR_cutoff = 0.75,
Plot_score_abs_cutoff = -0.1,
...
)

```

Arguments

datafile	the data files' path for the analysis, leave it as blank to enable a graphical user interface to select the data
projectfolder	optional, if NULL script will extract the path from datafile(s), and use the first workdir as project folder
threshold	specify the intensities threshold (0 to 1 in percentage) to report a identified molecule
ppm	the mz tolerance (in ppm) for peak integration
Digestion_site	Set the enzyme digestion specificity by one or more regex expressions or the name of a enzyme
missedCleavages	miss cleavage number allowed in this PMF search
Fastadatabase	the fasta database used in this pmf search, the file should be placed in the same folder with data files
adducts	the adducts list to be used for generating the PMF search candidates
Modifications	set the modifications

Substitute_AA	set the amino acid Substitutions
Decoy_search	enable (default) or disable the decoy search
Decoy_adducts	define the adduct list for decoy search. the decoy adducts could be "M+ACN+H", "M+IsoProp+H", "M+Rn".
Decoy_mode	select the decoy search mode between "isotope" (default), "element" and "adduct"
mzrange	define the mz range for the experiment, default is 700 to 4000 m/z.
IMS_analysis	Set "true" if you want to perform data pre-processing and proteomics search, set "false" if you want to bypass it
peptide_ID_filter	set the minimal count of peptides needed to identify a protein
Protein_feature_summary	"IMS_analysis" follow-up process that will collect all the identified peptide information and associate them with possible proteins
Peptide_feature_summary	"IMS_analysis" follow-up process that will summarize all datafiles identified peptides and generats a "peptide shortlist" in the result summary folder
plot_ion_image	"Peptide_feature_summarya" follow-up process that will plot every connponents in the "peptide shortlist". please use the cluster image grid to output the images.
parallel	the number of threads will be used in the PMF search, this option now only works for windows OS
spectra_segments_per_file	optimal number of distinctive regions in the tissue section, a virtual segmentation will be applied to the image files with this value. To have a better PMF result you may set a value that in the sweet point of sensitivety and false discovery rate (FDR).
Segmentation	set as "spatialKMeans" to enable a "spatialKMeans" Segmentation; set as "spatialShrunkenCentroids" to enable a "spatialShrunkenCentroids" Segmentation; If a region rank file was supplied, you can set this as "Virtual_segmentation" to perform a manual segmentation; Set it as "none" to bypass the segmentation.
preprocess	a list of params that define the IMS data pre-processing procedure
Smooth_range	"Segmentation" pixel smooth range
Virtual_segmentation_rankfile	specify a region rank file contains region information for manually region segmentation
Rotate_IMG	specify a configuration file to further change the rotation of the images
Region_feature_summary	"IMS_analysis" follow-up process that will summarize mz feature of all regions of all data files into the summary folder
use_previous_candidates	set as TRUE to reload the previously generated candidate list.
score_method	specify the peptide spectrum scoring method, "SQ RTP" is recommended.
plot_cluster_image_grid	set as "TRUE" to enable the protein cluster image function.
plot_cluster_image_overwrite	Set as true to generate the cluster images regardless the existance of previously file(s)

componentID_colname	Specify the component ID column in the result spreadsheet.
ClusterID_colname	Specify the cluster ID column in the result spreadsheet.
Protein_desc_of_interest	Specify a list of protein descriptions for cluster image plotting. Default setting will plot all reported proteins.
Protein_desc_of_exclusion	Specify a list of protein descriptions to be excluded from cluster image plotting.
plot_unique_component	Set as "TRUE" to plot only the unique components in the cluster image plotting.
FDR_cutoff	set the protein FDR cutoff threshold, default is 5 percent
plot_matching_score	enable the spectrum matching overlay plot
Component_plot_coloure	set as "mono" to use a pre-defined color scale to plot component images. Set as "as.cluster" to use the previously assigned mono color in the additive cluster binning process.
cluster_color_scale	Set as "blackwhite" to use only black and white color in the cluster image plotting. using "blackwhite" in cluster_color_scale will overwrite the components' color setting.
plot_layout	Set as "line" to plot cluster and component images for multiple data file or as "grid" to plot cluster images for single data file. In "grid" mode, Image's will be rendered into a grid with 5 columns.
export_Header_table	Set as "TRUE" to plot the header in the cluster image plotting. Header table includes the basic information of cluster and components.
export_footer_table	Set as "TRUE" to plot the footer in the cluster image plotting. Footer shows the protein coverage in the Proteomics mode.
attach_summary_cluster	Set as "TRUE" to attach an enlarged cluster image to the bottom of the cluster image.
remove_cluster_from_grid	Set as "TRUE" to remove the cluster image from the cluster image grid. it is recommended to set this same as the attach_summary_cluster.
cluster_rds_path	set as NULL if there is not preprocessed.rds available for a single file, script will load the raw data file which may reduce the signal intensities. For multiple samples, scripts will try to load the RDS file from each "ID" folder and merge the mz features via instrument resolution setting and output a combined RDS file to the project folder. For multiple files cluster images rendering user should set the attach_summary_cluster as False, and set remove_cluster_from_grid as true.

Value

None

Examples

```
imaging_identification(threshold=0.05, ppm=5,Digestion_site="[G]",
                        missedCleavages=0:1,Fastadatabase="murine_matrisome.fasta",
                        adducts=c("M+H", "M+NH4", "M+Na"),IMS_analysis=TRUE,
                        Protein_feature_summary=TRUE,plot_cluster_image=TRUE,
                        Peptide_feature_summary=TRUE,plot_ion_image=FALSE,
                        parallel=3,spectra_segments_per_file=5,Segmentation="spatialKMeans"
                        )
```

imaging_Spatial_Quant *imaging_Spatial_Quant*

Description

This is a spatial quantitation function for maldi imaging data set this function will read the candidate list file and generate quantification result

Usage

```
imaging_Spatial_Quant(
  datafile = tk_choose.files(filter = matrix(c("imzml file", ".imzML", "Text", ".txt",
    "All files", "*"), 3, 2, byrow = TRUE), caption =
    "Choose single or multiple file(s) for analysis"),
  threshold = 0,
  ppm = 2.5,
  Quant_list = "Metabolites of Interest.csv",
  adducts = c("M-H", "M+Cl"),
  cal.mz = F,
  mzlist_bypass = T,
  IMS_analysis = TRUE,
  Protein_feature_summary = T,
  plot_cluster_image = T,
  plot_style = "fleximaging",
  Peptide_feature_summary = T,
  plot_ion_image = FALSE,
  parallel = detectCores()/2,
  spectra_segments_per_file = 5,
  Smooth_range = 1,
  Segmentation = c("spatialKMeans", "spatialShrunkenCentroids", "Virtual_segmentation",
    "none"),
  Virtual_segmentation_rankfile = tk_choose.files(default =
    "Z:/George skyline results/maldiimaging/Maldi_imaging - Copy/radius_rank.csv",
    caption = "Choose Virtual segmentation rank info file"),
  Spectrum_feature_summary = T,
  Region_feature_summary = T,
  Region_feature_analysis = T,
  plot_each_metabolites = T,
  Cluster_level = "High",
  ClusterID_colname = "Name",
  Region_feature_analysis_bar_plot = T,
```

```

    norm_datafiles = T,
    norm_Type = "Median",
    mzrange = "auto-detect",
    BPPARAM = bpparam(),
    Rotate_IMG = NULL,
    ...
)

```

Arguments

datafile	specify the imzML data files
threshold	specify the intensities threshold (0 to 1 in percentage) to report a identified molecule
ppm	the mz tolerance (in ppm) for peak integration
Quant_list	the quantification candidate list, spatial quantification will go through every datafile and collect the ion intensities for each listed component
adducts	the adducts list to be used for generating the PMF search candidates
cal.mz	If set with "true", the function will recalculate the mz value according to the column named "formular" in the Quant_list and the specified adducts.
mzlist_bypass	Set "true" if you want to bypass the mzlist generating process
Protein_feature_summary	"IMS_analysis" follow-up process that will collect all the identified peptide information and associate them with possible proteins
plot_cluster_image	"Protein_feature_summary" follow-up process that will plot the protein cluster image
plot_ion_image	"Peptide_feature_summary" follow-up process that will plot every components in the "peptide shortlist"
parallel	the number of threads will be used in the PMF search, this option now only works for windows OS
spectra_segments_per_file	optimal number of distinctive regions in the imaging, a virtual segmentation will be applied to the image files with this value. To have a better PMF result you may set a value that in the sweet point of sensitivity and false discovery rate (FDR).
Smooth_range	"Segmentation" pixel smooth range
Segmentation	set as "spatialKMeans" to enable a "spatialKMeans" Segmentation; set as "spatialShrunkenCentroids" to enable a "spatialShrunkenCentroids" Segmentation; If a region rank file was supplied, you can set this as "Virtual_segmentation" to perform a manual segmentation; Set it as "none" to bypass the segmentation.
Virtual_segmentation_rankfile	specify a region rank file contains region information for manually region segmentation
Peptide_feature_summary	"IMS_analysis" follow-up process that will summarize all datafiles identified peptides and generates a "peptide shortlist" in the result summary folder

Value

None

Examples

```
imaging_Spatial_Quant(threshold=0.05, ppm=5,Digestion_site="[G]",
                        missedCleavages=0:1,Fastadatabase="murine_matrisome.fasta",
                        adducts=c("M+H", "M+NH4", "M+Na"),IMS_analysis=TRUE,
                        Protein_feature_summary=TRUE,plot_cluster_image=TRUE,
                        Peptide_feature_summary=TRUE,plot_ion_image=FALSE,
                        parallel=3,spectra_segments_per_file=5,spatialKMeans=TRUE
                        )
```

merge_atoms	<i>merge_atoms</i>
-------------	--------------------

Description

This is a function that merge the atoms to a molecule atom list .

Usage

```
merge_atoms(
  atoms,
  addelements,
  check_merge = F,
  mode = c("add", "ded"),
  multiplier = c(1, 1)
)
```

Arguments

atoms	the file name of candidate list
addelements	the folder that contains candidate list
check_merge	Check the final formula for any negative number of atom(s).
mode	add the elements to the formula by using "add", and deduct the elements by using "ded".
multiplier	Two number to define the atom and adding elements' proportion respectively.

Value

an atom list

Examples

```
merge_atoms(atoms=get_atoms("C7H13O8P"),addelements=list(H=1),check_merge=T,mode="ded",multiplier=c(1,2))
```

Meta_feature_list_fun *Meta_feature_list_fun*

Description

This is a function that prepare the candiate list for maldi imaging data qualitative or quantitative analysis. this function will read the candidate list file and generate mz for the adducts list defined in "adducts".

Usage

```
Meta_feature_list_fun(  
  database,  
  workdir = getwd(),  
  adducts = c("M-H", "M+Cl"),  
  cal.mz = TRUE,  
  bypass = FALSE,  
  BPPARAM = bpparam()  
)
```

Arguments

database	the file name of candidate list
workdir	the folder that contains candidate list
adducts	the adducts list to be used for generating the PMF search candidates
cal.mz	If set with "TRUE", the function will recalculate the mz value according to the column named "formula" in the database and the specified adducts.
BPPARAM	parallel processing parameter for BiocParallel
mzlist_bypass	Set "TRUE" if you want to bypass the mzlist generating process, the function will keep the mz and adduct as it is for the furture analysis. Be sure that the table contains "mz", "adduct" and "moleculeNames" as they are essential for later steps.

Value

a table of candiate list

Examples

```
Meta_feature_list_fun(database="lipid_candidates.csv",adducts=c("M-H","M+Cl"))
```

Peptide_modification	<i>Peptide_modification</i>
----------------------	-----------------------------

Description

This is a function that prepare the modification list for maldi imaging data qualitative or quantitative analysis. this function will load a pre-built Unimod database (unimod.df) into global environment.

Usage

```
Peptide_modification(retrieve_ID = NULL, mod_position = NULL, update_unimod = F)
```

Arguments

retrieve_ID	the file name of candidate list
mod_position	the folder that contains candidate list
update_unimod	the adducts list to be used for generating the PMF search candidates

Value

a table of modification list

Examples

```
Peptide_modification()
```

Protein_feature_list_fun	<i>Protein_feature_list_fun</i>
--------------------------	---------------------------------

Description

This is a function that prepare the candidate list for maldi imaging data qualitative or quantitative analysis. this function will read the fasta file and generate mz for the adducts list defined in "adducts".

Usage

```
Protein_feature_list_fun(
  workdir = getwd(),
  database,
  Digestion_site = "[G]",
  missedCleavages = 0:1,
  Multiple_mode = c("sequential", "parallel"),
  adducts = c("M+H", "M+NH4", "M+Na"),
  BPPARAM = bpparam(),
  Decoy_search = T,
  Decoy_mode = c("adducts", "elements", "isotope", "sequence"),
```

```

Decoy_adducts = c("M+He", "M+Ne", "M+Ar", "M+Kr", "M+Xe", "M+Rn"),
Substitute_AA = list(AA = c(NULL), AA_new_formula = c(NULL), Formula_with_water =
  c(NULL)),
mzrange = c(500, 4000),
output_candidatelist = T,
Modifications = list(fixed = NULL, fixmod_position = NULL, variable = NULL,
  varmod_position = NULL),
use_previous_candidates = F,
Protein_desc_of_exclusion = NULL,
Database_stats = F
)

```

Arguments

<code>workdir</code>	the folder that contains fasta file
<code>database</code>	the file name of fasta file
<code>Digestion_site</code>	Digestion pattern, this is the
<code>missedCleavages</code>	Define a number range for misscleaved peptides to be considered
<code>adducts</code>	the adducts list to be used for generating the PMF search candidates
<code>BPPARAM</code>	parallel processing parameter for BiocParallel

Value

a table of peptide candiate list

Examples

```
Protein_feature_list_fun(database="lipid candidates.csv",adducts=c("M+H","M+NH4","M+Na"))
```

Index

Cleavage_rules_fun, [2](#)
cluster_image_grid, [3](#)

get_atoms, [3](#)
grid.ftable, [4](#)

HiTMaP_GUI, [4](#)

imaging_identification, [5](#)
imaging_Spatial_Quant, [9](#)

merge_atoms, [11](#)
Meta_feature_list_fun, [12](#)

Peptide_modification, [13](#)
Protein_feature_list_fun, [13](#)